



Cranial Remolding Orthoses

CAD/CAM Technology

Becker Orthopedic has offered central fabrication services for **cranial remolding orthoses** for over six years using traditional plaster molds. In an effort to enhance this service, Becker is pleased to offer its customers the option of having these devices fabricated from electronic files utilizing the OMEGA® Tracer® software and OMEGA® Scanner.*



Orthotists can now quickly capture the shape of an infant's cranium using the hand-held OMEGA® scanner and no longer need to worry about the challenges associated with casting infants. Once the shape is captured electronically, the orthotist can then send the shape file (modified or unmodified) with a completed order form to Becker Orthopedic for fabrication.

 **OMEGA**
Tracer



**To inquire about OMEGA® Tracer® software and equipment, please call Ohio Willow Wood at (800) 848-4930 or visit the OWW website at owwco.com.*